

Product datasheet for **GA114528**

Human NIRF (UHRF2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UHRF2 CRISPRa kit - CRISPR gene activation of human ubiquitin like with PHD and ring finger domains 2
Symbol:	UHRF2
Components:	GA114528G1, NIRF gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114528G2, NIRF gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114528G3, NIRF gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96PU4</u>
Synonyms:	NIRF, URF2, RNF107, MGC33463, DKFZp434B0920, DKFZp686G0837, RP11-472F14.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)