

Product datasheet for **GA114534**

Human SLC5A11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC5A11 CRISPRa kit - CRISPR gene activation of human solute carrier family 5 member 11
Symbol:	SLC5A11
Components:	GA114534G1, SLC5A11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114534G2, SLC5A11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114534G3, SLC5A11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WWX8
Synonyms:	KST1; RKST1; SGLT6; SMIT2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)