

Product datasheet for **GA114535**

Human BAFF Receptor (TNFRSF13C) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF13C CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 13C
Symbol:	TNFRSF13C
UniProt ID:	Q96RJ3
Synonyms:	BAFF-R; BAFFR; BROMIX; CD268; CVID4; prolixin
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)