

## Product datasheet for **GA114536**

### Human Killer cell immunoglobulin like receptor 3DL3 (KIR3DL3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KIR3DL3 CRISPRa kit - CRISPR gene activation of human killer cell immunoglobulin like receptor, three Ig domains and long cytoplasmic tail 3
Symbol:	KIR3DL3
Synonyms:	CD158Z; KIR3DL7; KIR44; KIRC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)