

Product datasheet for **GA114546**

Human C3ORF24 (FANCD2OS) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FANCD2OS CRISPRa kit - CRISPR gene activation of human FANCD2 opposite strand
Symbol:	FANCD2OS
Components:	GA114546G1, C3ORF24 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114546G2, C3ORF24 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114546G3, C3ORF24 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96PS1
Synonyms:	C3orf24
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)