

Product datasheet for **GA114553**

Human C16orf42 (TSR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSR3 CRISPRa kit - CRISPR gene activation of human TSR3 ribosome maturation factor
Symbol:	TSR3
Components:	GA114553G1, C16orf42 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114553G2, C16orf42 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114553G3, C16orf42 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9UJK0
Synonyms:	C16orf42
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)