

Product datasheet for **GA114554**

Human ODAD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCDC151 CRISPRa kit – CRISPR gene activation of human coiled-coil domain containing 151
Symbol:	ODAD3
Components:	GA114554G1, ODAD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114554G2, ODAD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114554G3, ODAD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A5D8V7</u>
Synonyms:	CILD30
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

