

Product datasheet for **GA114563**

Human DNTTIP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNTTIP1 CRISPRa kit - CRISPR gene activation of human deoxynucleotidyltransferase terminal interacting protein 1
Symbol:	DNTTIP1
UniProt ID:	Q9H147
Synonyms:	C20orf167; dj447F3.4; Tdif1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)