

Product datasheet for **GA114566**

Human LRP15 (LRRC3B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC3B CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 3B
Symbol:	LRRC3B
Components:	GA114566G1, LRP15 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114566G2, LRP15 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114566G3, LRP15 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96PB8</u>
Synonyms:	LRP15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)