

Product datasheet for **GA114567**

Human KLHDC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KLHDC3 CRISPRa kit - CRISPR gene activation of human kelch domain containing 3
Symbol:	KLHDC3
Components:	GA114567G1, KLHDC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114567G2, KLHDC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114567G3, KLHDC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BQ90
Synonyms:	dj20C7.3; hPEAS; PEAS; RP1-20C7.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)