

Product datasheet for **GA114580**

Human **TLCD1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TLCD1 CRISPRa kit - CRISPR gene activation of human TLC domain containing 1
Symbol:	TLCD1
Components:	GA114580G1, TLCD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114580G2, TLCD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114580G3, TLCD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96CP7</u>
Synonyms:	TLC domain containing 1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)