

Product datasheet for **GA114585**

Human Pannexin 3 (PANX3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PANX3 CRISPRa kit - CRISPR gene activation of human pannexin 3
Symbol:	PANX3
UniProt ID:	<u>Q96QZ0</u>
Synonyms:	PX3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)