

Product datasheet for **GA114591**

Human TM4SF18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TM4SF18 CRISPRa kit - CRISPR gene activation of human transmembrane 4 L six family member 18
Symbol:	TM4SF18
Components:	GA114591G1, TM4SF18 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114591G2, TM4SF18 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114591G3, TM4SF18 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96CE8</u>
Synonyms:	L6D
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)