

Product datasheet for **GA114597**

Human CLNK activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CLNK CRISPRa kit - CRISPR gene activation of human cytokine dependent hematopoietic cell linker
Symbol:	CLNK
Components:	GA114597G1, CLNK gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114597G2, CLNK gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114597G3, CLNK gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z7G1</u>
Synonyms:	MIST
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)