

Product datasheet for **GA114598**

Human C1orf19 (TSEN15) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSEN15 CRISPRa kit - CRISPR gene activation of human tRNA splicing endonuclease subunit 15
Symbol:	TSEN15
UniProt ID:	Q8WW01
Synonyms:	C1orf19; sen15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)