

Product datasheet for **GA114616**

Human **WDR17** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	WDR17 CRISPRa kit - CRISPR gene activation of human WD repeat domain 17
Symbol:	WDR17
Components:	GA114616G1, WDR17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114616G2, WDR17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114616G3, WDR17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IZU2</u>
Synonyms:	FLJ26618
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

