

Product datasheet for **GA114618**

Human ACAP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ACAP3 CRISPRa kit - CRISPR gene activation of human ArfGAP with coiled-coil, ankyrin repeat and PH domains 3
Symbol:	ACAP3
UniProt ID:	Q96P50
Synonyms:	CENTB5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

