

Product datasheet for **GA114623**

Human CENTG3 (AGAP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AGAP3 CRISPRa kit - CRISPR gene activation of human ArfGAP with GTPase domain, ankyrin repeat and PH domain 3
Symbol:	AGAP3
UniProt ID:	Q96P47
Synonyms:	AGAP-3; CENTG3; cnt-g3; CRAG; MRIP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

