

Product datasheet for **GA114624**

Human TADA1L (TADA1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TADA1 CRISPRa kit - CRISPR gene activation of human transcriptional adaptor 1
Symbol:	TADA1
UniProt ID:	Q96BN2
Synonyms:	ADA1; hADA1; HFI1; STAF42; TADA1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)