

Product datasheet for **GA114629**

Human **SCGB3A2** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SCGB3A2 CRISPRa kit - CRISPR gene activation of human secretoglobin family 3A member 2
Symbol:	SCGB3A2
UniProt ID:	Q96PL1
Synonyms:	LU103; pnSP-1; PNSP1; UGRP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)