

Product datasheet for **GA114632**

Human WFIKKN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WFIKKN1 CRISPRa kit - CRISPR gene activation of human WAP, follistatin/kazal, immunoglobulin, kunitz and netrin domain containing 1
Symbol:	WFIKKN1
Components:	GA114632G1, WFIKKN1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114632G2, WFIKKN1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114632G3, WFIKKN1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96NZ8</u>
Synonyms:	C16orf12; RJD2; WFDC20A; WFIKKN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)