

## Product datasheet for **GA114633**

### Human RAB3IP activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	RAB3IP CRISPRa kit - CRISPR gene activation of human RAB3A interacting protein
Symbol:	RAB3IP
UniProt ID:	<u><a href="#">Q96QF0</a></u>
Synonyms:	RABIN3; RABIN8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)