


Product datasheet for **GA114642**

Human IHPK3 (IP6K3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IP6K3 CRISPRa kit - CRISPR gene activation of human inositol hexakisphosphate kinase 3
Symbol:	IP6K3
Components:	GA114642G1, IHPK3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114642G2, IHPK3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114642G3, IHPK3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96PC2</u>
Synonyms:	IHPK3; INSP6K3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)