

Product datasheet for **GA114645**

Human TAGAP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAGAP CRISPRa kit - CRISPR gene activation of human T cell activation RhoGTPase activating protein
Symbol:	TAGAP
UniProt ID:	Q8N103
Synonyms:	ARHGAP47; FKSG15; IDDM21; TAGAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)