

Product datasheet for **GA114646**

Human TMC1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMC1 CRISPRa kit - CRISPR gene activation of human transmembrane channel like 1
Symbol:	TMC1
UniProt ID:	Q8TDI8
Synonyms:	DFNA36; DFNB7; DFNB11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)