

Product datasheet for **GA114672**

Human **PIK3AP1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3AP1 CRISPRa kit - CRISPR gene activation of human phosphoinositide-3-kinase adaptor protein 1
Symbol:	PIK3AP1
UniProt ID:	Q6ZUJ8
Synonyms:	BCAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)