

Product datasheet for **GA114686**

Human C10orf78 (SFR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SFR1 CRISPRa kit - CRISPR gene activation of human SWI5 dependent homologous recombination repair protein 1
Symbol:	SFR1
Components:	GA114686G1, C10orf78 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114686G2, C10orf78 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114686G3, C10orf78 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q86XK3
Synonyms:	bA373N18.1; C10orf78; MEI5; MEIR5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)