

Product datasheet for **GA114688**

Human CLRN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CLRN3 CRISPRa kit - CRISPR gene activation of human clarin 3
Symbol:	CLRN3
Components:	GA114688G1, CLRN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114688G2, CLRN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114688G3, CLRN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NCR9</u>
Synonyms:	TMEM12; USH3AL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)