

Product datasheet for **GA114690**

Human PNLIPRP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PNLIPRP3 CRISPRa kit - CRISPR gene activation of human pancreatic lipase related protein 3
Symbol:	PNLIPRP3
UniProt ID:	<u>Q17RR3</u>
Synonyms:	OTTHUMP00000020563; pancreatic lipase-related protein 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)