

Product datasheet for **GA114709**

Human **SLC36A4** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | SLC36A4 CRISPRa kit - CRISPR gene activation of human solute carrier family 36 member 4 |
| Symbol: | SLC36A4 |
| UniProt ID: | Q6YBV0 |
| Synonyms: | PAT4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)