

## Product datasheet for **GA114737**

### Human LRIG3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LRIG3 CRISPRa kit - CRISPR gene activation of human leucine rich repeats and immunoglobulin like domains 3
Symbol:	LRIG3
UniProt ID:	<u><a href="#">Q6UXM1</a></u>
Synonyms:	LIG3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)