

Product datasheet for **GA114738**

Human TMEM132D activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM132D CRISPRa kit - CRISPR gene activation of human transmembrane protein 132D
Symbol:	TMEM132D
Components:	GA114738G1, TMEM132D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114738G2, TMEM132D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114738G3, TMEM132D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14C87</u>
Synonyms:	MOLT; PPP1R153
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)