

Product datasheet for **GA114750**

Human IKIP (IKBIP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IKBIP CRISPRa kit - CRISPR gene activation of human IKBKB interacting protein
Symbol:	IKBIP
Components:	GA114750G1, IKIP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114750G2, IKIP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114750G3, IKIP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q70UQ0</u>
Synonyms:	IKIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)