

Product datasheet for **GA114767**

Human PRR20A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRR20A CRISPRa kit - CRISPR gene activation of human proline rich 20A
Symbol:	PRR20A
Components:	GA114767G1, PRR20A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114767G2, PRR20A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114767G3, PRR20A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P86496</u>
Synonyms:	PRR20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)