

Product datasheet for **GA114786**

Human PPIL5 (LRR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRR1 CRISPRa kit - CRISPR gene activation of human leucine rich repeat protein 1
Symbol:	LRR1
Components:	GA114786G1, PPIL5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114786G2, PPIL5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114786G3, PPIL5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96L50
Synonyms:	4-1BBLRR; LRR-1; PPIL5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)