


Product datasheet for **GA114798**

Human **SLC24A4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC24A4 CRISPRa kit - CRISPR gene activation of human solute carrier family 24 member 4
Symbol:	SLC24A4
Components:	GA114798G1, SLC24A4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114798G2, SLC24A4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114798G3, SLC24A4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NFF2</u>
Synonyms:	AI2A5; NCKX4; SHEP6; SLC24A2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)