

Product datasheet for **GA114822**

Human DCUN1D3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DCUN1D3 CRISPRa kit - CRISPR gene activation of human defective in cullin neddylation 1 domain containing 3
Symbol:	DCUN1D3
Components:	GA114822G1, DCUN1D3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114822G2, DCUN1D3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114822G3, DCUN1D3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8IWE4
Synonyms:	44M2.4; SCCRO3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)