


Product datasheet for **GA114824**

Human CMTM3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CMTM3 CRISPRa kit - CRISPR gene activation of human CKLF like MARVEL transmembrane domain containing 3
Symbol:	CMTM3
Components:	GA114824G1, CMTM3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114824G2, CMTM3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114824G3, CMTM3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96MX0</u>
Synonyms:	BNAS2; chemokine-like factor super family 3; chemokine-like factor superfamily 3; CKLF-like MARVEL transmembrane domain containing 3; CKLFSF3; FLJ31762; MGC51956; OTTHUMP00000174790
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)