

Product datasheet for **GA114839**

Human ZNF720 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF720 CRISPRa kit - CRISPR gene activation of human zinc finger protein 720
Symbol:	ZNF720
Components:	GA114839G1, ZNF720 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114839G2, ZNF720 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114839G3, ZNF720 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7Z2F6
Synonyms:	MGC9313; MGC62019
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)