

Product datasheet for **GA114865**

Human WFIKKN2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WFIKKN2 CRISPRa kit – CRISPR gene activation of human WAP, follistatin/kazal, immunoglobulin, kunitz and netrin domain containing 2
Symbol:	WFIKKN2
Components:	GA114865G1, WFIKKN2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114865G2, WFIKKN2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114865G3, WFIKKN2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TEU8</u>
Synonyms:	GASP-1; hGASP-1; WFDC20B; WFIKKNRP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

