

Product datasheet for GA114876

Human SPNS2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPNS2 CRISPRa kit - CRISPR gene activation of human sphingolipid transporter 2
Symbol:	SPNS2
UniProt ID:	<u>Q8IVW8</u>
Synonyms:	spinster homolog 2; spinster homolog 2 (Drosophila)
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)