

Product datasheet for **GA114879**

Human WDR81 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR81 CRISPRa kit - CRISPR gene activation of human WD repeat domain 81
Symbol:	WDR81
Components:	GA114879G1, WDR81 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114879G2, WDR81 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114879G3, WDR81 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q562E7
Synonyms:	CAMRQ2; PPP1R166; SORF-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)