

Product datasheet for **GA114884**

Human keratin 40 (KRT40) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | KRT40 CRISPRa kit - CRISPR gene activation of human keratin 40 |
| Symbol: | KRT40 |
| UniProt ID: | Q6A162 |
| Synonyms: | CK-40; K40; KA36 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)