

Product datasheet for **GA114891**

Human **INO80C** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	INO80C CRISPRa kit - CRISPR gene activation of human INO80 complex subunit C
Symbol:	INO80C
Components:	GA114891G1, INO80C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114891G2, INO80C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114891G3, INO80C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PI98</u>
Synonyms:	C18orf37; hles6; IES6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

