

Product datasheet for **GA114896**

Human ZNF543 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF543 CRISPRa kit - CRISPR gene activation of human zinc finger protein 543
Symbol:	ZNF543
Components:	GA114896G1, ZNF543 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114896G2, ZNF543 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114896G3, ZNF543 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q08ER8
Synonyms:	DKFZp434H055; MGC119382; MGC119384
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)