

Product datasheet for **GA114910**

Human ZNF813 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF813 CRISPRa kit – CRISPR gene activation of human zinc finger protein 813
Symbol:	ZNF813
Components:	GA114910G1, ZNF813 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114910G2, ZNF813 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114910G3, ZNF813 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZN06</u>
Synonyms:	FLJ16542; FLJ34141
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

