

Product datasheet for **GA114911**

Human ZNF441 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF441 CRISPRa kit – CRISPR gene activation of human zinc finger protein 441
Symbol:	ZNF441
Components:	GA114911G1, ZNF441 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114911G2, ZNF441 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114911G3, ZNF441 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N8Z8</u>
Synonyms:	FLJ38637
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

