

Product datasheet for **GA114924**

Human ZNF787 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF787 CRISPRa kit - CRISPR gene activation of human zinc finger protein 787
Symbol:	ZNF787
Components:	GA114924G1, ZNF787 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114924G2, ZNF787 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114924G3, ZNF787 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6DD87</u>
Synonyms:	TIP20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)