

Product datasheet for **GA114926**

Human WDR88 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR88 CRISPRa kit - CRISPR gene activation of human WD repeat domain 88
Symbol:	WDR88
Components:	GA114926G1, WDR88 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114926G2, WDR88 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114926G3, WDR88 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZMY6</u>
Synonyms:	PQWD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

