

Product datasheet for **GA114927**

Human TMIGD2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMIGD2 CRISPRa kit - CRISPR gene activation of human transmembrane and immunoglobulin domain containing 2
Symbol:	TMIGD2
Components:	GA114927G1, TMIGD2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114927G2, TMIGD2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114927G3, TMIGD2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96BF3
Synonyms:	CD28H; IGPR-1; IGPR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)